

VOC Standard Solution - 33 components (EPA 8010B), 200 µg/mL of each component in Methanol (purge and trap)

Art. ID LM24-OP-4002165
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	alpha-Chlorotoluene	[100-44-7]	200	µg/mL		
	Bromobenzene	[108-86-1]	200	µg/mL		
	Bromoform (Tribromomethane)	[75-25-2]	200	µg/mL		
	Bromomethane	[74-83-9]	200	µg/mL		
	Carbon tetrachloride (Tetrachloromethane)	[56-23-5]	200	µg/mL		
	Chlorobenzene	[108-90-7]	200	µg/mL		
	Chloroethane	[75-00-3]	200	µg/mL		
	Chloroform (Trichloromethane)	[67-66-3]	200	µg/mL		
	Chloromethane	[74-87-3]	200	µg/mL		
	Chlorodibromomethane	[124-48-1]	200	µg/mL		
	Dibromomethane	[74-95-3]	200	µg/mL		
	1,2-Dichlorobenzene	[95-50-1]	200	µg/mL		
	1,3-Dichlorobenzene	[541-73-1]	200	µg/mL		
	1,4-Dichlorobenzene	[106-46-7]	200	µg/mL		
	Bromodichloromethane	[75-27-4]	200	µg/mL		
	Dichlorodifluoromethane	[75-71-8]	200	µg/mL		
	1,1-Dichloroethane	[75-34-3]	200	µg/mL		
	1,2-Dichloroethane	[107-06-2]	200	µg/mL		
	1,1-Dichloroethene	[75-35-4]	200	µg/mL		
	trans-1,2-Dichloroethene	[156-60-5]	200	µg/mL		
	1,2-Dichloropropane	[78-87-5]	200	µg/mL		
	cis-1,3-Dichloropropene	[10061-01-5]	200	µg/mL		
	trans-1,3-Dichloropropene	[10061-02-6]	200	µg/mL		
	Methylene chloride (Dichloromethane)	[75-09-2]	200	µg/mL		

chloromethane)				
1,1,1,2-Tetrachloroethane	[630-20-6]	200		µg/mL
1,1,2,2-Tetrachloroethane	[79-34-5]	200		µg/mL
Tetrachloroethene	[127-18-4]	200		µg/mL
1,1,1-Trichloroethane	[71-55-6]	200		µg/mL
1,1,2-Trichloroethane	[79-00-5]	200		µg/mL
Trichloroethene	[79-01-6]	200		µg/mL
Fluorotrichloromethane	[75-69-4]	200		µg/mL
1,2,3-Trichloropropane	[96-18-4]	200		µg/mL
Vinyl chloride	[75-01-4]	200		µg/mL