

ACs Internal Standard Solution - 2 components (EPA 8020A), 2000 µg/mL of each component in Methanol

Art. ID LM24-OP-4002176
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	alpha,alpha,alpha-Trifluorotoluene	[98-08-8]	2000	µg/mL		
	2-Bromofluorobenzene	[1072-85-1]	2000	µg/mL		