

Phenols Surrogate Standard Mixture (EPA 8041) - 2 components, 2000 µg/mL of each component in 2-Propanol

Art. ID LM24-OP-4002191
Unit 1 mL
Deliverydetails De Minimis UN1219 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Fluorophenol	[367-12-4]	2000	µg/mL		
	2,4,6-Tribromophenol	[118-79-6]	2000	µg/mL		