

PCBs Standard Solution - 19 components (EPA 8082), 100 µg/mL of each component in 2,2,4-Trimethylpentane (Isooctane)

Art. ID LM24-OP-4002212
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Chlorobiphenyl (PCB 1)	[2051-60-7]	100	µg/mL		
	2,3-Dichlorobiphenyl (P CB 5)	[16605-91-7]	100	µg/mL		
	2,2',5-Trichlorobipheny I (PCB 18)	[37680-65-2]	100	µg/mL		
	2,4',5-Trichlorobipheny I (PCB 31)	[16606-02-3]	100	µg/mL		
	2,2',3,5'-Tetrachlorobi phenyl (PCB 44)	[41464-39-5]	100	µg/mL		
	2,2',5,5'-Tetrachlorobi phenyl (PCB 52)	[35693-99-3]	100	µg/mL		
	2,3',4,4'-Tetrachlorobi phenyl (PCB 66)	[32598-10-0]	100	µg/mL		
	2,2',3,4,5'-Pentachloro biphenyl (PCB 87)	[38380-02-8]	100	µg/mL		
	2,2',4,5,5'-Pentachloro biphenyl (PCB 101)	[37680-73-2]	100	µg/mL		
	2,3,3',4',6-Pentachloro biphenyl (PCB 110)	[38380-03-9]	100	µg/mL		
	2,2',3,4,4',5'-Hexachlo robiphenyl (PCB 138)	[35065-28-2]	100	µg/mL		
	2,2',3,4,5,5'-Hexachlor obiphenyl (PCB 141)	[52712-04-6]	100	µg/mL		
	2,2',3,5,5',6-Hexachlor obiphenyl (PCB 151)	[52663-63-5]	100	µg/mL		
	2,2',4,4',5,5'-Hexachlo robiphenyl (PCB 153)	[35065-27-1]	100	µg/mL		
	2,2',3,3',4,4',5-Heptac	[35065-30-6]	100	µg/mL		

chlorobiphenyl (PCB 170)				
2,2',3,4,4',5,5'-Heptac	[35065-29-3]	100		µg/mL
chlorobiphenyl (PCB 180)				
2,2',3,4,4',5,6-Heptac	[52663-69-1]	100		µg/mL
chlorobiphenyl (PCB 183)				
2,2',3,4',5,5',6-Heptac	[52663-68-0]	100		µg/mL
chlorobiphenyl (PCB 187)				
2,2',3,3',4,4',5,5',6-N	[40186-72-9]	100		µg/mL
nonachlorobiphenyl (PCB 206)				