

ACs Surrogate Standard Mixture - 6 components, 4000 µg/mL of each component in Methylene chloride (Dichloromethane)

Art. ID LM24-OP-4002265
Unit 1 mL
Deliverydetails De Minimis UN1593 III E1

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Fluorobiphenyl	[321-60-8]	4000	µg/mL		
	2-Fluorophenol	[367-12-4]	4000	µg/mL		
	Nitrobenzene D5	[4165-60-0]	4000	µg/mL		
	Phenol D5	[4165-62-2]	4000	µg/mL		
	p-Terphenyl D14	[1718-51-0]	4000	µg/mL		
	2,4,6-Tribromophenol	[118-79-6]	4000	µg/mL		