

**PAH Internal Standard Mixture - 6 components (EPA 8275A), 4000 µg/mL of each component in Methylene chloride (Dichloromethane)**

Art. ID LM24-OP-4002268  
Unit 1 mL  
Deliverydetails De Minimis UN1593 III E1

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene D10	[15067-26-2]	4000	µg/mL		
	Chrysene D12	[1719-03-5]	4000	µg/mL		
	1,4-Dichlorobenzene D4	[3855-82-1]	4000	µg/mL		
	Naphthalene D8	[1146-65-2]	4000	µg/mL		
	Perylene D12	[1520-96-3]	4000	µg/mL		
	Phenanthrene D10	[1517-22-2]	4000	µg/mL		