

ACs Standard Solution - 6 components (ISO 11423-2,1997) 50 ml/L of each component in Acetone

Art. ID LM24-OP-4002344
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	50	mL/L		
	Toluene	[108-88-3]	50	mL/L		
	o-Xylene	[95-47-6]	50	mL/L		
	m-Xylene	[108-38-3]	50	mL/L		
	p-Xylene	[106-42-3]	50	mL/L		
	Ethylbenzene	[100-41-4]	50	mL/L		