

Internal and injection Standard 1 component (EN 15308:2008)
2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl (PCB 209) 400 µg/mL in n-Hexane

Art. ID LM24-OP-4002371
Unit 1 mL
Deliverydetails De Minimis UN1208 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl (P CB 209)	[2051-24-3]	400	µg/mL		