

Calibration Mixture without PCB 170 (acc. method 2) (EN 12766), 13 components, 10 µg/mL of each component in 2,2,4-Trimethylpentane (Isooctane)

Art. ID LM24-OP-4002372
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',5-Trichlorobiphenyl (PCB 18)	[37680-65-2]	10	µg/mL		
	2,4,4'-Trichlorobiphenyl (PCB 28)	[7012-37-5]	10	µg/mL		
	2,4',5-Trichlorobiphenyl (PCB 31)	[16606-02-3]	10	µg/mL		
	2,2',3,5'-Tetrachlorobiphenyl (PCB 44)	[41464-39-5]	10	µg/mL		
	2,2',5,5'-Tetrachlorobiphenyl (PCB 52)	[35693-99-3]	10	µg/mL		
	2,2',4,5,5'-Pentachlorobiphenyl (PCB 101)	[37680-73-2]	10	µg/mL		
	2,3',4,4',5-Pentachlorobiphenyl (PCB 118)	[31508-00-6]	10	µg/mL		
	2,2',3,4,4',5'-Hexachlorobiphenyl (PCB 138)	[35065-28-2]	10	µg/mL		
	2,2',3,4',5',6-Hexachlorobiphenyl (PCB 149)	[38380-04-0]	10	µg/mL		
	2,2',4,4',5,5'-Hexachlorobiphenyl (PCB 153)	[35065-27-1]	10	µg/mL		
	2,2',3,4,4',5,5'-Heptachlorobiphenyl (PCB 180)	[35065-29-3]	10	µg/mL		
	2,2',3,3',4,4',5,5'-Octachlorobiphenyl (PCB 194)	[35694-08-7]	10	µg/mL		
	2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl (PCB 209)	[2051-24-3]	10	µg/mL		