

**EN 14039 - Retention time window Standard Solution - 2 components, 30 mg/L of each component in n-Heptane**

Art. ID LM24-OP-4002429  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1206 II E2

Description

Validity: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Decane (C10)	[124-18-5]	30	mg/L		
	n-Tetracontane (C40)	[4181-95-7]	30	mg/L		