

CEN/TS 16621, PAH4 standard Solution - of 4 components, 10 µg/mL of each component in Acetonitrile

Art. ID LM24-OP-4002452
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

Description

Validity: 12 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------|------------|---------------------|-------|--------|--------|
| | Benzo(a)pyrene | [50-32-8] | 10 | µg/mL | | |
| | Benzo(a)anthracene | [56-55-3] | 10 | µg/mL | | |
| | Benzo(b)fluoranthene | [205-99-2] | 10 | µg/mL | | |
| | Chrysene | [218-01-9] | 10 | µg/mL | | |