

Surrogate Mixture (Method 8260) - 3 components in Methanol (purge & trap) 2500 µg/mL each

Art. ID LM24-OP-4002482
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dibromofluoromethane	[1868-53-7]	2500	µg/mL		
	Toluene D8	[2037-26-5]	2500	µg/mL		
	4-Bromofluorobenzene	[460-00-4]	2500	µg/mL		