

Multiresidue Pesticide Standard (Mix 3) - 25 components, 100 µg/mL each in Toluene / Acetonitrile (99/1)

Art. ID LM24-OP-4002498
Unit 1 mL
Deliverydetails De Minimis UN1294 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benfluralin	[1861-40-1]	100	µg/mL		
	Biphenyl	[92-52-4]	100	µg/mL		
	Chlorothalonil	[1897-45-6]	100	µg/mL		
	Dichlofluanid	[1085-98-9]	100	µg/mL		
	Dicloran	[99-30-9]	100	µg/mL		
	3,4-Dichloroaniline	[95-76-1]	100	µg/mL		
	Dichlobenil	[1194-65-6]	100	µg/mL		
	Diphenylamine	[122-39-4]	100	µg/mL		
	Ethalfuralin	[55283-68-6]	100	µg/mL		
	Fluchloralin	[33245-39-5]	100	µg/mL		
	Isopropalin	[33820-53-0]	100	µg/mL		
	Nitralin	[4726-14-1]	100	µg/mL		
	Nitrofen	[1836-75-5]	100	µg/mL		
	Oxyfluorfen	[42874-03-3]	100	µg/mL		
	Pendimethalin	[40487-42-1]	100	µg/mL		
	Pentachloroaniline	[527-20-8]	100	µg/mL		
	Pentachlorobenzonitrile	[20925-85-3]	100	µg/mL		
	Pentachloronitrobenzene	[82-68-8]	100	µg/mL		
	Prodiamine	[29091-21-2]	100	µg/mL		
	Profluralin	[26399-36-0]	100	µg/mL		
	2,3,5,6-Tetrachloroaniline	[3481-20-7]	100	µg/mL		
	Tecnazene	[117-18-0]	100	µg/mL		
	cis-1,2,3,6-Tetrahydrophthalimide (THPI)	[1469-48-3]	100	µg/mL		
	Tolyfluanid	[731-27-1]	100	µg/mL		
	Trifluralin	[1582-09-8]	100	µg/mL		