

Standard mixture - 8 components in Iso-octane

Art. ID LM24-OP-4002507
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|--------------|---------------------|------|--------|--------|
| | 2,4,4'-Trichlorobiphenyl (PCB 28) | [7012-37-5] | 10 | mg/L | | |
| | 2,2',5,5'-Tetrachlorobiphenyl (PCB 52) | [35693-99-3] | 10 | mg/L | | |
| | 2,4',5-Trichlorobiphenyl (PCB 31) | [16606-02-3] | 10 | mg/L | | |
| | 2,2',4,5,5'-Pentachlorobiphenyl (PCB 101) | [37680-73-2] | 10 | mg/L | | |
| | 2,3',4,4',5-Pentachlorobiphenyl (PCB 118) | [31508-00-6] | 10 | mg/L | | |
| | 2,2',3,4,4',5'-Hexachlorobiphenyl (PCB 138) | [35065-28-2] | 10 | mg/L | | |
| | 2,2',4,4',5,5'-Hexachlorobiphenyl (PCB 153) | [35065-27-1] | 10 | mg/L | | |
| | 2,2',3,4,4',5,5'-Heptachlorobiphenyl (PCB 180) | [35065-29-3] | 10 | mg/L | | |