

**Semi-Volatiles GC/MS Tuning Standard - 4 components, 1000 µg/mL of each component in Methylene chloride (Dichloromethane)**

Art. ID LM24-OP-4002560  
Unit 1 mL  
Deliverydetails De Minimis UN1593 III E1

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bis(pentafluorophenyl)phosphine	[5074-71-5]	1000	µg/mL		
	4,4'-Benzidine	[92-87-5]	1000	µg/mL		
	Pentachlorophenol	[87-86-5]	1000	µg/mL		
	4,4'-DDT	[50-29-3]	1000	µg/mL		