

Base/Neutrals Surrogate Standard Mixture - 3 components, 1000 µg/mL of each component in Methylene chloride (Dichloromethane)

Art. ID LM24-OP-4002567
Unit 1 mL
Deliverydetails De Minimis UN1593 III E1

Description

Validity: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | 2-Fluorobiphenyl | [321-60-8] | 1000 | µg/mL | | |
| | Nitrobenzene D5 | [4165-60-0] | 1000 | µg/mL | | |
| | p-Terphenyl D14 | [1718-51-0] | 1000 | µg/mL | | |