

PAH Internal Surrogate Standard Fortification Solution - (EPA 525.2) - 3 components, 500 µg/mL of each component in Acetone

Art. ID LM24-OP-4002574
Unit 1 mL
Deliverydetails De Minimis UN1090 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene D10	[15067-26-2]	500	µg/mL		
	Chrysene D12	[1719-03-5]	500	µg/mL		
	Phenanthrene D10	[1517-22-2]	500	µg/mL		