

Method DM 471 Standard Mixture - 1 - 8 components, 1000 µg/mL of each component in Methanol

Art. ID LM24-OP-4002647
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	1000	µg/mL		
	Ethylbenzene	[100-41-4]	1000	µg/mL		
	Styrene	[100-42-5]	1000	µg/mL		
	Methyl-tert-butylether (MTBE)	[1634-04-4]	1000	µg/mL		
	Toluene	[108-88-3]	1000	µg/mL		
	o-Xylene	[95-47-6]	1000	µg/mL		
	m-Xylene	[108-38-3]	1000	µg/mL		
	p-Xylene	[106-42-3]	1000	µg/mL		