

Method DM 471 Standard Mixture 3 - 7 components, 10 µg/mL of each component in Methanol

Art. ID LM24-OP-4002649
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|------------|---------------------|-------|--------|--------|
| | Aniline | [62-53-3] | 10 | µg/mL | | |
| | Diphenylamine | [122-39-4] | 10 | µg/mL | | |
| | 2-Amino-1-methylbenzene (o-Toluidine) | [95-53-4] | 10 | µg/mL | | |
| | 2-Anisidine | [90-04-0] | 10 | µg/mL | | |
| | 3-Anisidine | [536-90-3] | 10 | µg/mL | | |
| | 4-Anisidine | [104-94-9] | 10 | µg/mL | | |
| | p-Toluidine | [106-49-0] | 10 | µg/mL | | |