

**Method DM 471 Standard Mixture 4 - 8 components, 500 µg/mL of each component in Methanol**

Art. ID LM24-OP-4002650  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Pentane	[109-66-0]	500	µg/mL		
	Hexane (Solvent blank)	[110-54-3]	500	µg/mL		
	n-Heptane (C7)	[142-82-5]	500	µg/mL		
	n-Octane (C8)	[111-65-9]	500	µg/mL		
	n-Nonane (C9)	[111-84-2]	500	µg/mL		
	n-Decane (C10)	[124-18-5]	500	µg/mL		
	n-Undecane (C11)	[1120-21-4]	500	µg/mL		
	n-Dodecane (C12)	[112-40-3]	500	µg/mL		