

**Method DM 471 Standard Mixture 5 - 7 components, 100 µg/mL of each component in Methanol**

Art. ID LM24-OP-4002651  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aniline	[62-53-3]	100	µg/mL		
	Diphenylamine	[122-39-4]	100	µg/mL		
	2-Amino-1-methylbenzene (o-Toluidine)	[95-53-4]	100	µg/mL		
	2-Anisidine	[90-04-0]	100	µg/mL		
	3-Anisidine	[536-90-3]	100	µg/mL		
	4-Anisidine	[104-94-9]	100	µg/mL		
	p-Toluidine	[106-49-0]	100	µg/mL		