

**Method DM 471 Standard Mixture - 10 - 7 components, 100 µg/mL of each component in Methanol**

Art. ID LM24-OP-4002654  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,1,2-Trichloroethane	[79-00-5]	100	µg/mL		
	Trichloroethene	[79-01-6]	100	µg/mL		
	1,2,3-Trichloropropane	[96-18-4]	100	µg/mL		
	1,1,2,2-Tetrachloroethane	[79-34-5]	100	µg/mL		
	Tetrachloroethene	[127-18-4]	100	µg/mL		
	Hexachloro-1,3-butadiene	[87-68-3]	100	µg/mL		
	1,1,1-Trichloroethane	[71-55-6]	100	µg/mL		