

Deuterated PAH Mixture - 3 components, 200 µg/mL of each component in Methylene chloride (Dichloromethane)

Art. ID LM24-OP-4002667
Unit 1 mL
Deliverydetails De Minimis UN1593 III E1

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthylene D8	[93951-97-4]	200	µg/mL		
	Chrysene D12	[1719-03-5]	200	µg/mL		
	Indeno(1,2,3-c,d)pyrene D12	[203578-33-0]	200	µg/mL		