

**Deuterated PAH Mixture - 3 components, 200 µg/mL of each component in Methylene chloride (Dichloromethane)**

Art. ID LM24-OP-4002667  
Unit 1 mL  
Deliverydetails De Minimis UN1593 III E1

Description

Validity: 24 months

| Text/Information | Analyte/Parameter              | CAS number    | Concentration/Value | Unit  | Method | Source |
|------------------|--------------------------------|---------------|---------------------|-------|--------|--------|
|                  | Acenaphthylene D8              | [93951-97-4]  | 200                 | µg/mL |        |        |
|                  | Chrysene D12                   | [1719-03-5]   | 200                 | µg/mL |        |        |
|                  | Indeno(1,2,3-c,d)pyrene<br>D12 | [203578-33-0] | 200                 | µg/mL |        |        |