

Deuterated PAH Mixture - 7 components, 200 µg/mL of each component in Methylene chloride (Dichloromethane)

Art. ID LM24-OP-4002668
Unit 1 mL
Deliverydetails De Minimis UN1593 III E1

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene D10	[15067-26-2]	200	µg/mL		
	Benzo(a)anthracene D12	[1718-53-2]	200	µg/mL		
	Benzo(a)pyrene D12	[63466-71-7]	200	µg/mL		
	Dibenz(a,h)anthracene D 14	[13250-98-1]	200	µg/mL		
	Dibenzo(a,i)pyrene-D14	[158776-07-9]	200	µg/mL		
	Fluoranthene D10	[93951-69-0]	200	µg/mL		
	Phenanthrene D10	[1517-22-2]	200	µg/mL		