

Deuterated PAH Mixture - 7 components, 200 µg/mL of each component in Methylene chloride (Dichloromethane)

Art. ID LM24-OP-4002668
Unit 1 mL
Deliverydetails De Minimis UN1593 III E1

Description

Validity: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------|---------------|---------------------|-------|--------|--------|
| | Acenaphthene D10 | [15067-26-2] | 200 | µg/mL | | |
| | Benzo(a)anthracene D12 | [1718-53-2] | 200 | µg/mL | | |
| | Benzo(a)pyrene D12 | [63466-71-7] | 200 | µg/mL | | |
| | Dibenz(a,h)anthracene D 14 | [13250-98-1] | 200 | µg/mL | | |
| | Dibenzo(a,i)pyrene-D14 | [158776-07-9] | 200 | µg/mL | | |
| | Fluoranthene D10 | [93951-69-0] | 200 | µg/mL | | |
| | Phenanthrene D10 | [1517-22-2] | 200 | µg/mL | | |