

PAH Calibration Mixture - 16 components in Acetonitrile

Art. ID LM24-OP-4002670
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene	[83-32-9]	20	µg/mL		
	Acenaphthylene	[208-96-8]	15	µg/mL		
	Anthracene	[120-12-7]	0,8	µg/mL		
	Benzo(a)anthracene	[56-55-3]	4	µg/mL		
	Benzo(a)pyrene	[50-32-8]	5	µg/mL		
	Benzo(b)fluoranthene	[205-99-2]	4	µg/mL		
	Benzo(g,h,i)perylene	[191-24-2]	3,5	µg/mL		
	Benzo(k)fluoranthene	[207-08-9]	4,5	µg/mL		
	Chrysene	[218-01-9]	3,5	µg/mL		
	Dibenz(a,h)anthracene	[53-70-3]	3,5	µg/mL		
	Fluoranthene	[206-44-0]	8	µg/mL		
	Fluorene	[86-73-7]	5	µg/mL		
	Indeno(1,2,3-cd)pyrene	[193-39-5]	4,5	µg/mL		
	Naphthalene	[91-20-3]	20	µg/mL		
	Phenanthrene	[85-01-8]	3,5	µg/mL		
	Pyrene	[129-00-0]	8,5	µg/mL		