

EPA 8010 Halogenated VOC Mix - 14 components, 2000 µg/mL of each component in Methanol

Art. ID LM24-OP-4002685
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bromoform (Tribromomethane)	[75-25-2]	2000	µg/mL		
	Carbon tetrachloride (Tetrachloromethane)	[56-23-5]	2000	µg/mL		
	Chloroform (Trichloromethane)	[67-66-3]	2000	µg/mL		
	Dibromomethane	[74-95-3]	2000	µg/mL		
	1,1-Dichloroethane	[75-34-3]	2000	µg/mL		
	1,2-Dichloroethane	[107-06-2]	2000	µg/mL		
	1,2-Dichloropropane	[78-87-5]	2000	µg/mL		
	1,1,1,2-Tetrachloroethane	[630-20-6]	2000	µg/mL		
	1,1,2,2-Tetrachloroethane	[79-34-5]	2000	µg/mL		
	Tetrachloroethene	[127-18-4]	2000	µg/mL		
	1,1,1-Trichloroethane	[71-55-6]	2000	µg/mL		
	1,1,2-Trichloroethane	[79-00-5]	2000	µg/mL		
	Trichloroethene	[79-01-6]	2000	µg/mL		
	1,2,3-Trichloropropane	[96-18-4]	2000	µg/mL		