

ICP MIX 12 components in 2 %HNO3

Art. ID LM24-OP-4002700
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	3000	mg/L		
	Magnesium (Mg)	[7439-95-4]	450	mg/L		
	Potassium (K)	[7440-09-7]	1200	mg/L		
	Sodium (Na)	[7440-23-5]	3000	mg/L		
	Sulfate (SO42-)		3000	mg/L		
	Phosphorus (P)	[7723-14-0]	234,8	mg/L		
	Iron (Fe)	[7439-89-6]	30	mg/L		
	Manganese (Mn)	[7439-96-5]	30	mg/L		
	Copper (Cu)	[7440-50-8]	15	mg/L		
	Zinc (Zn)	[7440-66-6]	15	mg/L		
	Boron (B)	[7440-42-8]	30	mg/L		
	Molybdenum (Mo)	[7439-98-7]	6	mg/L		