

**Tuning Solution for ICP/MS 9 components 10 mg/L each in HNO3 2% (Standard equivalent to Agilent/Varian)**

Art. ID LM24-OP-4002702  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	10	mg/L		
	Beryllium (Be)	[7440-41-7]	10	mg/L		
	Cerium (Ce)	[7440-45-1]	10	mg/L		
	Cobalt (Co)	[7440-48-4]	10	mg/L		
	Indium (In)	[7440-74-6]	10	mg/L		
	Magnesium (Mg)	[7439-95-4]	10	mg/L		
	Lead (Pb)	[7439-92-1]	10	mg/L		
	Thorium (Th)	[7440-29-1]	10	mg/L		
	Thallium (Tl)	[7440-28-0]	10	mg/L		