

EN 16995: 2017 - Internal MOSH/MOAH Standard Solution 1- 9 components in Toluene

Art. ID LM24-OP-4002703
Unit 1 mL
Deliverydetails De Minimis UN1294 II E2

Description

Shelf-Life: 12 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Undecane (C11)	[1120-21-4]	300	µg/ml		
	n-Tridecane (C13)	[629-50-5]	150	µg/ml		
	1,3,5-Tri-tert-butylbenzene	[1460-02-2]	300	µg/ml		
	Dicyclohexyl	[92-51-3]	300	µg/ml		
	1-Methylnaphthalene	[90-12-0]	300	µg/ml		
	2-Methylnaphthalene	[91-57-6]	300	µg/ml		
	n-Pentylbenzene	[538-68-1]	300	µg/ml		
	Perylene	[198-55-0]	600	µg/ml		
	alpha,alpha,alpha-20R-c holestane	[481-21-0]	600	µg/ml		