

**ICP Multi-Element Standard - EPA 200.7 - Mixed Calibration Standard 5 - 6 components in 5 %HNO<sub>3</sub>**

Art. ID LM24-OP-4002715

Unit 100 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	50	mg/L		
	Nickel (Ni)	[7440-02-0]	20	mg/L		
	Beryllium (Be)	[7440-41-7]	10	mg/L		