

## **Quality Control - 23 components (Standard Equivalent to Jobin Yvon)**

Art. ID LM24-OP-4002936  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5] | 1000                | mg/L |        |        |
|                  | Barium (Ba)       | [7440-39-3] | 1000                | mg/L |        |        |
|                  | Bismuth (Bi)      | [7440-69-9] | 1000                | mg/L |        |        |
|                  | Boron (B)         | [7440-42-8] | 1000                | mg/L |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 1000                | mg/L |        |        |
|                  | Calcium (Ca)      | [7440-70-2] | 1000                | mg/L |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 1000                | mg/L |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 1000                | mg/L |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 1000                | mg/L |        |        |
|                  | Gallium (Ga)      | [7440-55-3] | 1000                | mg/L |        |        |
|                  | Indium (In)       | [7440-74-6] | 1000                | mg/L |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 1000                | mg/L |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 1000                | mg/L |        |        |
|                  | Lithium (Li)      | [7439-93-2] | 1000                | mg/L |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 1000                | mg/L |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 1000                | mg/L |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 1000                | mg/L |        |        |
|                  | Potassium (K)     | [7440-09-7] | 1000                | mg/L |        |        |
|                  | Silver (Ag)       | [7440-22-4] | 1000                | mg/L |        |        |
|                  | Sodium (Na)       | [7440-23-5] | 1000                | mg/L |        |        |
|                  | Strontium (Sr)    | [7440-24-6] | 1000                | mg/L |        |        |
|                  | Thallium (Tl)     | [7440-28-0] | 1000                | mg/L |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 1000                | mg/L |        |        |