

Standard for Semi-quantitative Method - 7 components (Standard Equivalent to Jobin Yvon)

Art. ID LM24-OP-4002938
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium (K)	[7440-09-7]	1000	mg/L		
	Silver (Ag)	[7440-22-4]	100	mg/L		
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Barium (Ba)	[7440-39-3]	100	mg/L		
	Boron (B)	[7440-42-8]	100	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		
	Silicon (Si)	[7440-21-3]	50	mg/L		