

Standard for the Determination of traces of 5 heavy metals (Standard Equivalent to Jobin Yvon)

Art. ID LM24-OP-4002942
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	50	mg/L		
	Selenium (Se)	[7782-49-2]	50	mg/L		
	Lead (Pb)	[7439-92-1]	30	mg/L		