

Standard Quality Control for ICP Ultima family - 5 components (Standard Equivalent to Jobin Yvon)

Art. ID LM24-OP-4002946
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium (K)	[7440-09-7]	1500	mg/L		
	Lead (Pb)	[7439-92-1]	1000	mg/L		
	Aluminium (Al)	[7429-90-5]	500	mg/L		
	Magnesium (Mg)	[7439-95-4]	500	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		