

**Standard Quality Control for Testing ICP Activa Family - 5 components (Standard Equivalent to Jobin Yvon)**

Art. ID LM24-OP-4002949  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		