

**Bismuth (Bi) 10000 mg/L in 5-10% HNO<sub>3</sub>, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4002955  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	10000	mg/L		Bi, 99.999 %