

## **Copper (Cu) 10000 mg/L in 2-5% HNO<sub>3</sub>, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4002960  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

### Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	10000	mg/L		Cu, 99.999 %