

Potassium (K) 10000 mg/L in 2-5% HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4002963
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------------------------------|
| | Potassium (K) | [7440-09-7] | 10000 | mg/L | | KNO ₃ , 99.9 99% |