

Manganese (Mn) 10000 mg/L in 2-5% HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4002966
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Manganese (Mn)	[7439-96-5]	10000	mg/L		Mn(NO ₃) ₂ , 99.999%