

Arsenic (As) 1000 mg/L in diluted HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4002976
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|-----------------|
| | Arsenic (As) | [7440-38-2] | 1000 | mg/L | | As, 99.999 % |