

Beryllium (Be) 1000 mg/L in diluted HNO3 for ICP/MS

Art. ID LM24-OP-4002981
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Beryllium (Be) | [7440-41-7] | 1000 | mg/L | | |