

Bismuth (Bi) 1000 mg/L in 2% HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4002982
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	1000	mg/L		Bi, 99.999 %