

## **Iron (Fe) 1000 mg/L in 2-5% HNO<sub>3</sub>, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4002988  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

### Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iron (Fe)	[7439-89-6]	1000	mg/L		Fe(NO <sub>3</sub> ) <sub>3</sub> , 99.9995%