

Germanium (Ge) 1000 mg/L in 2-5% HNO3/HF tr, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4002990
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Germanium (Ge)	[7440-56-4]	1000	mg/L		Ge, 99.999 %