

**Potassium (K) 1000 mg/L in 2-5% HNO<sub>3</sub>, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4002994  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium (K)	[7440-09-7]	1000	mg/L		KNO <sub>3</sub> , 99.9 99%