

Sodium (Na) 1000 mg/L in diluted HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003000
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium (Na)	[7440-23-5]	1000	mg/L		NaNO ₃ , 99.999%