

## **Phosphorus (P) 1000 mg/L in H<sub>2</sub>O, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4003002  
Unit 100 mL  
Deliverydetails No Dangerous Good /not restricted

### Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphorus (P)	[7723-14-0]	1000	mg/L		NH <sub>4</sub> H <sub>2</sub> PO <sub>4</sub> , 99.999%