

**Lead (Pb) 1000 mg/L in 2-5% HNO<sub>3</sub>, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4003003  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lead (Pb)	[7439-92-1]	1000	mg/L		Pb(NO <sub>3</sub> ) <sub>2</sub> , 99.9995%